

Thermo Scientific Orion

AquaMate Spectrophotometer External Device Interfacing and Data Transfer Guide

The following instructions describe how to interface a Thermo Scientific[™] Orion[™] AquaMate[™] spectrophotometer with an external device, transfer the data and save the data as a Microsoft[®] Excel[®], Word[®] or Notepad file. The operator must set the AquaMate spectrophotometer settings to allow data to be sent to an external device, install the AquaMate USB driver on the computer, and install a computer program for a data collection program such as Realterm or HyperTerminal.

AquaMate Spectrophotometer Setup for Computer Interfacing



- 1. Power on the Orion AquaMate spectrophotometer.
- 2. Press the *utility* key to acess the instrument utility menu.



3. Press the ▲ or ▼ key to highlight the **Printer** option and press the *enter* key to open the printer menu.

Utility	11:55am 23Jun14
Software Revision: Instrument Serial Number: Gringer Date/Time Setup Standby Baseline Expiration (hr:m Baseline Expiration (hr:m	1.004 2W2R091104 0991 11:55am 23Jun14 30 min 30 ff
Beeper Language % of lamp life used Screen Contrast Printout Contrast	On English 0.55%
Stored Tests Directory USB Drive Files Calculator <u>Press ↑ or ↓ to select pa</u>	urameter.

4. Press the ▲ or ▼ key to highlight the **PC** option and press the *enter* key to set the PC option to **On**.

Utility	11:56am 23Jun14
Software Revisi Instrument Seri	on: 1.004 al Number: 2W2R091104
Date/Time Setu Date/Time Setu Standby Baseline Expir Beeper Language X of lamp life Screen Contras	Internal Off PC Off External Off Press f or J to select Press ENTER to change Press ESC to exit
Printout Contr Stored Tests Di USB Drive Files Calculator	rectory
Press ↑ or ↓ to	select parameter.

5. Press the *esc* key until the instrument returns to the main measurement screen.



USB Computer Cable Interfacing and USB Driver Installation



USB communication port

- 1. Connect the USB computer cable to the USB communication port on the back of the AquaMate spectrophotometer.
- 2. Insert the Orion AquaMate Documentation DVD into the computer CD/DVD drive.
- 3. Connect the USB computer cable to a USB port on the computer.
- 4. Open the computer Device Manager.
 - a. Click the Start icon.



b. Right-click the **Computer** button.



c. Click the Properties button.



d. Click the Device Manager button.

000		to Warmen and the second	
Control Panel +	All Control Panel Items + System	 • • • • Secrit Control Parel 	
(in Edit View Tools theip			
Cantrol Fond Hone Cantrol Manager Cantrol Manager Cantrol Manager Cantrol Manager System protection Cantrol system settings	View basic information Windows edition Windows / Protesponi Copyright © 2009 Microw rights reserved. Service Pack I	about your computer)
	System		
	Ratings	Wedgest Expiriment Index	
	Processor	Intel(R) Core(TN) IS CPU M 526 @ 240GHz 240 GHz	
	Installed memory (RAM):	4.00 GB (3.24 GD usable)	
	Sustain turner	12-bit Operating System	
	where the		
Ser and Performance beformation and Local	Pen and Touch	No Pen or Touch Input is available for this Display	

5. Double-click the **Spectrophotometer** button on the Device Manager window.



6. Click the **Update Driver** button on the Spectrophotometer Properties window.

ectrophot	ometer Propert	ties	
Seneral D	river Details		
10	Spectrophotome	ter	
	Device type:	Other devices	
	Manufacturer	Unknown	
	Location	Port_#0003.Hub_#0004	
Device s	tatus		
The dri	vers for this devic	ce are not installed. (Code 28)	i.
There is	no driver select	ted for the device information set or element.	
To find	a driver for this d	evice, click Update Driver.	
			*
		Update Driver_	_
		have a second	_
		OK	Cance

7. Click the **Browse my computer for driver** software button.



8. Click the Browse button.



- 9. Select the appropriate USB driver folder from the Orion AquaMate Documentation DVD in the computer CD/DVD drive and click the OK button.
 - a. Select the **Vista32** folder when using a computer with 32-bit Microsoft 7 or Vista operating system.
 - b. Select the **Vista64** folder when using a computer with 64-bit Microsoft 7 or Vista operating system.
 - c. Select the **XP32** folder when using a computer with Microsoft XP operating system.

Browse For Folder	×
Select the folder that contains drivers for you	r hardware.
4 🎼 Computer	
> 🥾 OSDisk (C:)	
4 🦨 DVD RW Drive (D:) Thermo	
D Thermo Scientific	
Uista32	
Uista64	=
📜 XP32	-
Vieta22	
Eolder:	
ОК	Cancel

10. Click the Next button.



11. Wait while the USB driver software is located and installed.

😋 🛓 liptin Directofere - Spectrykolowie) 🚡 Shidate Diver Schware - Spectrophytometer	100
Searching DWistal2 for software_		Installing driver software	_
	Gener		

12. Approve any security messages by clicking the **Install** button.

Windows Security	×
Would you like to install this device softwa	ire?
Name: Thermo Scientific, Inc. Ports (COM & L	PT
Aways trust software from "Thermo Fisher Scientific Inc.".	install Dog't Install
You should only install driver software from publish device software is safe to install?	ers you trust. How can I decide which

13. Once the installation of the USB driver software is complete, click the **Close** button.



14. Record the COM port number of the Spectrophotometer USB connection in the Device Manager window.

E Device Manager	- • ×
Elle Action View Help	
今今(四)(日)(日四)(お)(日夜あ	
) 🥁 Batteries	
Bluetooth Radios	
🗈 🕷 Computer	
E ControlVault Device	
Disk drives	
Display adapters	
》 凸 DVD/CD-ROM drives	
» Q資 Human Interface Devices	
IDE ATA/ATAPI controllers	-
IEEE 1394 Bus host controllers	1
Keyboards	
If the and other pointing devices	
) 🖨 Moderns	
Monitors	
Network adapters	
Ports (COM & EPT)	
- TOP ECP Printer Port (LPT1)	
Spectrophotometer (COM9)	
Processors	
> go Smart card readers	

Example Data Collection Program

The AquaMate communcation settings are as follows:

- Bits per second: 19200
- Data bits: 8
- Parity: None
- Stop bits: 1
- Flow control: None
- Select a computer program for data collection and make sure the program, including device drivers and plug-ins, is fully installed. Realterm data collection software can be downloaded from their website and is used in the following instructions as an example.
- 2. Open the data collection software program by clicking the Realterm icon on the computer desktop.



3. Click the **Port** tab on the RealTerm window. Confirm that the RealTerm communication settings match the appropriate AquaMate communcation settings.

	a second starts				
0.4			an lana	Less Links	Maad Ersamel 1
isplay Por	Cepture	Pana Send Echo Port	ec 50%	1204sc 3]	Clear Freeze
isplay Pot	Cepture Eort [3	Pins Send Eche Port	12C 202	tCMsc h jn ⊈ Charge	Clear Freeze
Isplay Port	Ceprune Port 3 Deta Bits Biblio	Pine Send Educ Port	I2C 202	DAnge	Clear Freeze
Hiptey Por Jeud 57600 Pgnly ■ None C Odd C Even	Cepure Port 3 Deta Bits @ Bits C 7 bits	Pess Send Echo Port	Dpan Spy Software Roy Receive	DOMise 5 In (Control Xon Charge 7 Xof Char 17 Xof Char 19	Clear Freeze 2 Consect PAD (7) DOP (7) CTS (8) DCD (7)
Hipley Por Pank Fank F None C Odd C Even C Merk	Cepure Deta Bits © 8 bits C 7 bits C 6 bits	Pms Sand Echo Port	IEC 202 Open Sty Software Roy Receive Transmit	DOMise h]n] Charge ⁷⁷ Costrol Xon Char, 17 Xof Char, 19 Wesockis -	Electric 1 Electric 1 000 (2) 1 700 (2) 1 700 (2) 1 700 (2) 1 700 (2) 1 700 (2) 1 700 (2) 1 700 (2) 1 700 (2) 1 700 (2) 1 700 (2) 1 700 (2) 1 1 1 700 (2) 1 1 1 700 (2) 1 1 1 700 (2) 1 1 1 700 (2) 1 1 1 700 (2) 1 1 1 700 (2) 1 1 1 700 (2) 1 1 1 700 (2) 1 1 1 700 (2) 1 1 1 700 (2) 1 1 1 700 (2) 1 1 1 700 (2) 1 1 1 700
Party Port Party P	Cepure Port 3 Dece Bits © 8 bits C 7 bits C 5 bits	Pins Sand Echo Port	RC 202 Open Spy Satware Ro F Receive	Costal Costal Xoo Clex 17 Xof Cher 19 Wesock is Pape C Pape C Pape	Clear Freeze 1

4. Click the **Port** dropdown menu and select the COM port number of the Spectrophotometer USB connection.

	(C) (C) X
vaphy PM Dapwe FVs Since (cheffer) 20 2001 2034c) 445 Imp Since Since Since Cheffer 20 2001 2034c) 445 Party Since Since 2011	yn Clear Freeze -)

5. Click the **Baud** dropdown menu and select **19200**.

RealTerri: Serial Capture Program 2.0.070	100 (E) X
Corpuy Per Corpus Pra Stad (State) 20 (201 (2034c) Mc)	(n) Clear Freezel :
State State (State) Pro (State) P	

6. Click the **Change** button to save the port settings.



7. Click the **Display** tab on the RealTerm window. Select the **Ascii** option under the **Display As** field.



- 8. When data from the AquaMate spectrophotometer is ready to be sent to the computer, first test the connection between the spectrophotometer and data collection program.
 - a. Press the *print* key on the AquaMate spectrophotometer. The data from the AquaMate spectrophotometer should appear on the program display.

Statistics Rato Print Noto Eave Pata	- say to Prog	011 #/ 	001007
B(ar Pactor) C Broult Bolts Abs 108 S18.4nn	Research for	-11, 21, 24, 27, 24, 04 -11, 24, 24, 24, 24 -11, 24, 24, 24, 24 -14, 27, 24, 24, 04 -14, 27, 24, 24, 04 -14, 24, 24, 24, 24 -14, 24, 24, 24 -14, 24, 24, 24 -14, 24 -14	
1 0,000 2 0,000 3 0,000	8,838 9,839 9,839 9,839	66+ 181: 191 ⁻	
Display As Display As Without Fr C Ass C Participated	laptice Fins 	Send Edus Port EC EC E ECMac M	c <u>in Clear Freeze 1</u> 9akr

- 9. Once the connection has been validated, select the type, name and location of the file the data will be saved as.
 - a. Click the **Capture** tab on the RealTerm window.



b. Click the ... button next to the File field and set the file type, name and location. For example, use Microsoft Excel (.xls) file type, name the file AQ8000-Data.xls and save the file to the computer desktop.

Seve As						×
Save js	E Desktop		•	****		
Recet Places	My Documents	Computer				
	File game. Save as gype	AQ80006-Distraji ete			•	- <u>S</u> avo Carcel

10. Click the **Start: Overwrite** button on the Capture tab of the RealTerm window. The program is now ready to receive data from the AquaMate spectrophotometer and export the data to the selected file type, name and location.

RealTerm: Serial Captur	e Program 2.0.0.70					X
						-
		m.	_		100	
play Port Capture	ins Send Ech	o Pot I2C	120-2	In Clear	Freeze	7
pley Port Capture	ins Send Ech	io Poit I2C	120-2	\n Clear	Freeze Stetus	?] ;
play Port Capture P apture: Stort Overwrite	Ins Send Ech	io Pot I2C Stog Capture	120-2 End Afe	\n Clear	Status	? F
play Port Capture i phone Start Overwrite 1	ins Send Ech Inst Append 5	to Pot I2C	I I I I I I I I I I I I I I I I I I I	\n Clear	Freeze Status Discon	? nect
pley Port Capture 1 phre Start Overwrite 1 IC C\Usersfricki hartwei	ins Send Ech Ins Append S Desktop/A08000 Ce	io Pot 12C Stog Capture staxi 💌 🔔	End Ate	\n Clear	Status Status Discon PodD (2 TXD (3	*
play Port Caprure phree Start Overwrite IC (Users/picki hartwei Capture as Harx	ins Send Ech Int Append S DesktoplA08000-Ce	io Pot I2C Stog Capture Maxi 👤 💶	End After Back Development Ca	∖n Clear	Freeze	? = = = = = = = = = = = = = = = = = = =
oley Port Cepure phre Stat Overwrite IC C:Usersfricki hartweil Capture as Hox	ins Send Ech Int Append S Desktop/A08000-Ce	to Pot I2C Stog Capture Ita xi • - Trive Stamo • None	LEC2	\n Clear cture ₽ Detenter • come	Freeze	*
oley Port Capture s phate Start Overwrite Ie (C) Usersfrick i fortwei Capture as Hex	Ins Send Edit	to Port I2C Stog Capture Istaul •	IZO2	\n Clear cture ≠ Detenter comm	Freeze	*
olay Port Capture p phare Start Overwrite 3 Re CiUsers/nicki hortwel Capture as Hex	Ins Send Ech Itati Aspend S Desktop(A03000 Ce	to Pot 2C Stog Capture nts.xl ▼ Time Stants Nons Unar Unar	IZO2	\n Clear corror _ corror _ Corror _ Corror corror corror corror	Freeze	* ned*
splay Port Capture r aptive Start Overwite I ColUsers/incki.fortwell Copture as Hex	Ins Send Edu tert Append S Desktop/A03000 Ce	io Por 2C Stog Capture Istual Time Stamo Unav Unav Unav Unav	End Alter Dred Cs Dred Cs	(n) Clear cture 2 Detroite communication	Freeze Stetus PAD(2 TCD(2 DCD(1 DCD(1 DSR(6 Ring(9 BREAV	* * * * * * * * * * * * * * * * * * *

- 11. Press the *print* key on the AquaMate spectrophotometer to send the data from the instrument to the data collection program.
- 12. Once all data has been sent from the instrument to the computer and data collection program, click the **Stop Capture** button on the Capture tab of the RealTerm window.
- 13. Open the data file, for example AQ8000-Data.xls, and format the data as needed.

AquaMate Spectrophotometer Setup for Printer Interfacing

The AquaMate printer compatibility requirements are as follows:

- Printer must be PCL5 or PCL6
- Printer must be a printer only, not a multi-function printer
- For this example, a PCL5 HP Laser Jet 1300 printer is used
- HP 2055D and Office jet 6000 printers have been used in the past
- 1. Power on the Orion AquaMate spectrophotometer.



- 2. Press the *utility* key to acess the instrument utility menu.
- 3. Press the ▲ or ▼ key to highlight the **Printer** option and press the *enter* key to open the printer menu.



4. Press the ▲ or ▼ key to highlight the **External** option and press the *enter* key to set the External printer option to **On**.

Software Revisi	n: 1.0
Instrument Seri	I NUMBER: ZWZRU911
Date/Time Setu Standby Baseline Expir Beeper Language % of lamp life Screen Contras Printout Contr	Internal Off PC Off External Off Press f or 1 to select Press ENTER to change Press ESC to exit
Stored Tests Di USB Drive Files Calculator	ectory

5. Press the *esc* key until the instrument returns to the main measurement screen.



USB Printer Cable Interfacing and Printing Data

1. Connect the USB cable to the USB printer port on the back of the AquaMate spectrophotometer.



- 2. Connect the USB printer cable to the USB port on the printer.
- 3. When data from the AquaMate spectrophotometer is ready to be sent to the printer, press the *print* key on the AquaMate spectrophotometer.

thermoscientific.com/water

© 2014 Thermo Fisher Scientific Inc. All rights reserved. Microsoft, Excel and Word are registered trademarks of Microsoft Corporation. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries.

Water Analysis Instruments

North America Toll Free: 1-800-225-1480 Tel: 1-978-232-6000 info.water@thermo.com Netherlands Tel: (31) 020-4936270 info.water.uk@thermo.com

China Tel: (86) 21-68654588 wai.asia@thermofisher.com India Tel: (91) 22-4157-8800 wai.asia@thermofisher.com

Singapore Tel: (65) 6778-6876 wai.asia@thermofisher.com Japan Tel: (81) 045-453-9175 wai.asia@thermofisher.com

Australia Tel: (613) 9757-4300 in Australia (1300) 735-295 InfoWaterAU@thermofisher.com

