

November 18, 2014

Dear Valued Customer,

This letter is being provided to serve as notification of change related to the configuration of the Self Contained Biological Indicators (SCBIs) for use in steam sterilization processes your company has purchased from NAMSA. The product codes associated with the change include SCS-05 and SCS-06.


The product configuration will be changed in terms of the SCBI vial label. The vial label currently contains a ISO Class 1 Chemical Indicator (CI) which transitions from yellow to dark brown and consists of a laminated paper. The vial label CI will be changed to an un-laminated paper label with Class 1 CI that transitions from blue to pink. The new label will be removable after exposure or incubation.

The changes are being made to improve the label's CI stability related to light sensitivity.

The change will be phased in as the current inventory levels of the existing configuration are depleted. We anticipate the timeframe for transition to be 90 days.

If you have any questions or concerns regarding this change, please do not hesitate to contact our Customer Service Group at 800.860.1888 or productorders@namsa.com, Jennifer Rogers, Inside Sales at 419.662.4421 or jrogers@namsa.com or myself using the details below.

Best regards,



General Manager
T: 419.662.4488
E: jwheeler@namsa.com