



Spare Part List to Steroclave® Autoclave: Parts for Models 1915X, 1925X, 1941X and 25X

Part #	Description
	Bottom only (specify size and with or without parts)
	Cover only (specify size and with or without parts)
54	Clamp bolt (630000)
55	Pin for clamp bolt (630010)
64	Bakelite wing nut (630020)
65	Control valve (640000)
72S	Geared steam gauge (640610)
73	Lens for steam gauge (Replacement) (670050)
75	Excess pressure release plug (Replaced by No. 2050CS)
76	Bakelite top handle (630070)
77	Bakelite top handle screw (630100)
78	Retaining bayonet clamp (630090)
79	Retaining bayonet clamp screw (630100)
111	Rack for inside aluminum container for 1915X, 1925X, 25X (630110)
151	Rack for inside aluminum container for 1941X (630120)
405	Side handle complete for 1941X (630210)
1010	Overpressure plug for sterilizer – red color (630130)
2050CS	Excess pressure relief valve (640180)
2151	Stainless steel support stand (used inside 25X) (640020)
2153	On-Off toggle switch for 25X (640020)
2154A 120V	Pilot light for 25X-120 Volt models (640190)
2154B 240V	Pilot light for 25X-240 Volt models (640200)
2155-15	Air exhaust tube for 1915X (640210)
2155-25	Air exhaust tube for 1925X, 25X (640220)
2155-41	Air exhaust tube for 1941X (640230)
2156	Aluminum container for 25X (640030)
2157	Detachable power supply cord for 25X (for models made before 1980) (640040)
2157A 120V	Wiring harness, grounded 3-wire power supply cord for 25X-120V (640240)
2157B 240V	Wiring harness, grounded 3-wire power supply cord for 25X- 240V (640250)
2157 SRB	Strain relief bushing for 25X (640260)
*2158AK	Heating element with adaptor kit (for units made prior to 1984) (640270)
*2158AS 120V	Heating element for 25X- 120V..(1050 Watt)..(includes fiber gaskets) (640280)
*2158BS 240V	Heating element for 25X- 240V..(1050 Watt)..(includes fiber gaskets) (640290)
2158FG	Fiber gasket for #2158 heating element (two needed) (640310)
*2159	Thermostat for 25X (includes fiber gasket) (640050)
2159FG	Fiber gasket for 2159 thermostat (one needed) (640310)
2160	Bakelite heat control knob for thermostat for 25X (640060)
2161	Cast aluminum control box for 25X (Z002161.M)
2162	Aluminum container for 1925X (640070)
2163	Aluminum container for 1915X (640080)

2164	Aluminum container for 1941X (640090)
2180	Support base for 25X (will elevate unit 1 3/8" above table surface) (Z002180.S)
*6054	Thermometer, dual scale, C scale: 10-150, F scale: 50-300 (640140)

* We recommend that the Thermostat, Heating element and Thermometer be replaced at the factory.
NO GUARANTEE FOR RESPONSIBILITY FOR THE PROPER FUNCTIONING OF THESE PARTS
CAN BE ASSUMED BY THE COMPANY IF THEY ARE NOT INSTALLED AT THE FACTORY.