

GARVER MFG. INC.

Centrifuge Parts List:

Part #	Description
0001C	Large Bowl
0002C	Bowl Lid 4"-9" Glassware
0003C	Bowl Lid 2 ½"-6" Glassware
0004C	Baseplate Large
0005C	Small Bowl
0006C	Bowl lid with hinges
0007C	Baseplate small
0008C	4 Place 250ml rotor head
0009C	Double cup hanger
0010C	Single cup hanger
0011C	Cup
0012C	Single cup w/ hanger
0013C	Double cup w/ hanger
0014C	Bottle carrier small
0015C	Bottle carrier large
0016C	Brake disc & stem
0017C	Brake handle 7 ¾"
0018C	Brake handle 10"
0019C	Flat brass brake handle
0020C	Rotor Standard/Model #/SN
0021C	Rotor Super/Model #SN
0022C	Nut – Motor Top
0023C	Rotor Shaft
0024C	Stator Standard (Model #/SN)
0025C	Stator Super (Model #/SN)
0026C	Nut Motor bottom
0027C	Motor housing – Top
0028C	Nut for carrier & visiscope
0029C	Bearing Top of Motor
0030C	Bearing Bottom of Motor
0031C	Ball Bearing Set Screw
0032C	Ball Bottom of Motor
0033C	Motor housing bottom
0034C	Visiscope (Model #/SN)
0035C	Visiscope face
0036C	Visiscope Shaft
0037C	Visiscope Bushing Kit
0038C	Thermometer Holder
0039C	Thermometer Red Celsius
0040C	Pilot Light Assembly Small
0041C	Pilot Light Assembly Large

0042C	Control Knob
0043C	Toggle Switch
0044C	Ohmite Tap Switch
0045C	Thermostat
0046C	VS Control Box
0047C	VS Transformer 13537-UP
0048C	VS Transformer & light 1-13536
0049C	Heater 500W
0050C	Heater Pin (Set of 2)
0051C	Lava Bushing (Set of two)
0052C	Cord & Plug for heater
0053C	Bolt for heater clamp
0054C	Heater Clamp 500W
0055C	Heater 750W
0056C	Switch Box Triangle
0057C	Switch Box Small
0058C	Handle for lid
0059C	Bolts/bowl & mtr (Part #0001C)
0060C	Bolts/bowl & mtr (Part #0005C)
0061C	Bolt & bushing for thermostat
0062C	Bolt for attaching lid to bowl
0063C	Bolts & Nuts for heater guard
0064C	Bolts for bowl lid hinge
0065C	Hinge
0066C	Pilot light bulb
0067C	Bolt & Nut for Mtr Housing
0068C	Motor Super (Model #/SN)
0069C	Motor Standard (Model #/SN)
0070C	Cord w/ 3 prong plug
0071C	Cork Gasket
0072C	Rubber Plug
0073C	Milk Solid Adaptors
0074C	Gerber Adaptors
0075C	Rubber pads milk solids
0076C	Heater Terminal Insulator
0077C	2 ½' Heater wire
0078C	ON/OFF Plate
0079C	Serial # Plate
0080C	VS Serial # Plate
0081C	Pocket Tach
0082C	Babcock Video & Script
0083C	Mojonnier Video & Script
0084C	Warning Lid Sticker
0085C	Rubber Feet (set of 3)
0086C	Centrifuge Oil