

IFM Proficiency Testing Programs: Legionella					
Order Code	Program	Sample Prep	Available Test	Enroll By Date (mm/dd/yyyy)	Dispatch Date (mm/dd/yyyy)
LG1	Cooling Tower Water	<b>Sample Count:</b> 3 <b>Format:</b> Freeze-dried vial <b>Preparation:</b> Reconstitute to obtain 100mL	<b>Count (cfu/mL):</b> Plate count 21°C, Plate count 37°C Legionella pneumophila (SG1, SG2-14), Total Legionella pneumophila, Legionella spp, Total Legionella	09/15/2025	10/6/2025
LG2	Cooling Tower Water			02/10/2025	03/3/2025
LG2s	Cooling Tower Water (detection <i>by filtration methods</i> )	Sample Count: 3 Format: Freeze-dried vial Preparation: Reconstitute to obtain 2500mL	<b>Count (cfu/1000 mL):</b> Legionella pneumophila, Legionella species, Total Legionella, species and serogroup identification of targets	06/16/2025	07/7/2025
LG3	Cooling Tower Water	<b>Sample Count:</b> 3 <b>Format:</b> Freeze-dried vial <b>Preparation:</b> Reconstitute to obtain 100mL	<b>Count (cfu/mL):</b> Plate count 21°C, Plate count 37°C Legionella pneumophila (SG1, SG2-14), Total Legionella pneumophila, Legionella spp, Total Legionella	06/16/2025	07/7/2025