

Lab products

The **heart** of every lab

Up to **25%**
off



Save on **Nalgene Rapid-Flow** filtration products

Boost your productivity and save up to 25% on Thermo Scientific™ Nalgene™ Rapid-Flow™ Sterile Disposable Filter Units with PES membranes and Thermo Scientific™ Nalgene™ Rapid-Flow™ Sterile Disposable Bottle Top Filters with PES membranes. Get the support, performance and precision you need for a range of applications.

- Devices feature the fastest flow rate and lowest clogging of comparable alternatives
- 0.1 µm PES available for mycoplasma protection
- 0.2 µm PES devices are stem cell and mouse embryonic assay (MEA) tested
- Non-pyrogenic and endotoxin-free—Rapid-Flow filter units meet the most rigorous testing standards for vacuum filtration units

Save on products at the core of your workflows

Order using the catalog numbers below and use promo code **FILT2525** at checkout to take advantage of this offer.

Ordering information

Nalgene Rapid-Flow Sterile Disposable Filter Units with PES

Membrane diameter	Pore size	Volume	Quantity	Cat. No.
75 mm	0.1 µm	500 mL	12/case	566-0010
90 mm	0.1 µm	1,000 mL	12/case	567-0010
50 mm	0.2 µm	50 mL	12/case	564-0020
50 mm	0.2 µm	150 mL	12/case	565-0020
50 mm	0.2 µm	250 mL	12/case	568-0020
75 mm	0.2 µm	500 mL	12/case	566-0020
90 mm	0.2 µm	500 mL	12/Case	569-0020
90 mm	0.2 µm	1,000 mL	12/case	567-0020
50 mm	0.45 µm	150 mL	12/case	165-0045
50 mm	0.45 µm	250 mL	12/case	168-0045
75 mm	0.45 µm	500 mL	12/case	166-0045
90 mm	0.45 µm	500 mL	12/case	169-0045
90 mm	0.45 µm	1,000 mL	12/case	167-0045

Nalgene Rapid-Flow Sterile Disposable Bottle Top Filters with PES membranes

Pore size	Volume	Fits bottle neck size	Quantity	Cat. No.
0.2 µm	150 mL	45 mm	12/case	596-4520
0.2 µm	500 mL	33 mm	12/case	595-3320
0.2 µm	500 mL	45 mm	12/case	595-4520
0.2 µm	1,000 mL	45 mm	12/case	597-4520
0.45 µm	150 mL	45 mm	12/case	296-4545

 Contact your Thermo Fisher Scientific sales representative to learn more.

Discount valid on the purchase of list-price items only through thermofisher.com orders. No additional notice of applied discounts will be given after checkout. All promotions and discounts will be applied to eligible Thermo Fisher Scientific cart items only. Only one promo code can be applied per order.

This offer is valid January 1, 2025 through June 30, 2025.

Customer acknowledges that this offer may include a discount or other price reduction that must be properly and accurately accounted for and reported by customer in accordance with all federal and state laws, including without limitation the federal anti-kickback law (42 U.S.C. § 1320a-7b(b)(3)(A)) and regulations thereunder (42 C.F.R. §1001.952(h)).

For Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for customers' specific use or application. ©2024 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. STUDIO-24-198-3373