



Parker Balston Compressed Air Microbial Test Unit (CAMTU) Upgrade Program

The Balston CAMTU (Compressed Air Microbial Test Unit) is a lightweight, portable device used to verify the microbiological quality of compressed air in food plants and other GMP-sensitive environments. CAMTU supports cleanliness programs aligned with GMPs and SQF expectations by enabling rapid, on-site sampling where compressed air contacts product or product-contact surfaces. The Balston CAMTU kit is supplied as a complete field-use kit for compressed air microbial monitoring, including the CAMTU sampling device and accessories needed to connect to a compressed air source and perform sampling using agar plates.

To improve customer satisfaction, the CAMTU sampling housing has been updated to accept commercially available, ready-to-use 90mm agar plates. Users will now be able to source commercially prepared agar plates for use with the CAMTU kit. The updated kit is available as part # C01-0175, and the CAMTU housing only is available as part # C01-0176. The updated kit will not include agar plates.

Parker Hannifin
Industrial Gas Filtration and Generation Division
242 Neck Road | Haverhill, MA 01835 | +1.800.343.4048
parker.com/balstonfilters



Upgrade Program Offer:

To help customers with this product transition, we are offering a complimentary upgrade to new part # C01-0176. This program will be open until **8/31/2026**. The following criteria must be met for upgrade eligibility:

- The Balston CAMTU housing or kit shipped from Parker after **July 1, 2024**.
- Part numbers C01-0142 and C01-0136 (or variations of C01-0136) qualify - please use multiple forms for claiming different part numbers
- The information on this form must be provided complete and accurate

CAMTU UPGRADE FORM FOR ELIGIBLE PURCHASES

Part Number		Shipping Information:	
		Contact Name	
Quantity		Contact Phone	
Point of Purchase (Name of Authorized Parker Balston Distributor or Representative)		Business Name	
Original Purchase Information Needed from Authorized Parker Balston Distributor or Rep:		Shipping Address	
Purchase Order # (Must list all if multiple)			
Parker Sales Order # (Must list all if multiple)			

Parker Hannifin has final discretion to review all information provided and determine qualifying criteria for an upgrade. Filling out this form does not guarantee a complimentary CAMTU upgrade.